




High-fidelity bovine vulva model for teaching Day-1 competence to veterinary students

Susanna N. Taylora, Ling Yan Cheunga, Eloi R.G. Guarnieria



a Jockey Club College of Veterinary Medicine and Life Sciences, City University of Hong Kong, Hong Kong S.A.R.

Background

Buhner suturing is considered an essential skill for newly graduated farm animal veterinarians (1). Inspiration for this model was drawn from a presentation on an equine Caslick suture model (2).

Objective

To develop a high-fidelity bovine vulva model meeting the following key criteria:

- correct size and orientation
- realistic sensation of entering the subcutaneous layer
- A replica of a bovine vulva region was modelled freehand out of air-dry clay and sealed with acrylic primer.
- A reusable **mold** was created using Mold Star 15#.
- A **silicone** model was cast in layers, with each layer allowed to semicure:

a	Vaginal mucosa	Thickened% Ecoflex 0030#, blood colour*, 2mm, internal vertical surfaces
b	Skin	Ecoflex 0030#, brown colour*, 4mm, reinforced with strips of powermesh
C	Subcutaneous tissue	Exoflex Gel#, 10mm
d	Muscles	Ecoflex 0030#, blood colour*, 15mm

Ecoflex gel#, 5mm

Materials: # Smooth-On, Macungie, USA; % THI-VEXTM, Smooth-On, Macungie, USA; * SilcPigTM, Smooth-On, Macungie, USA

How to use this model?

Equipment

Bovine vulva model

Adhesive layer

- Buhner needle 12"
- Scalpel blade #15
- Umbilical tape 1/4"

Teaching material

Skill sheet

Instructions

- 1.Prepare equipment and teaching materials
- 2. Attach the model to a smooth vertical surface
- 3.Closely supervise students when using blades and needles, and specifically when exiting the skin with the Buhner needle
- 4.Be aware that excessively tightening the Buhner suture may cause the model to detach from the vertical surface
- 5.If required, repair the model (i.e., skin incisions) using additionally applied EcoflexTM 00-30A

Team's feedback

- We believe we successfully simulated the sensation of subcutaneous needle advancement.
- Our mold has produced four distinct models, demonstrating its replicability.
- While the gel layer temporarily secures the model to a smooth surface, it may weaken over time, and when force is applied with the needle.

References

1.Miller, R.B. et al., 2004. Practitioner-defined competencies required of new veterinary graduates in food animal practice. J. Vet. Med. Educ. 31, 347–365. https://doi.org/10.3138/jvme.31.4.347

2.Collingsworth, B., Hunt, J., 2022. Equine Caslick's model.

Dr. Eloi Guarnieri Clinical Assistant Professor (Ruminant Health Management)





